



Andersen EOGas® Sterilization bags require a heavy duty sealer to ensure good, consistent results.

The Andersen AN1091 heat sealer is a reliable, professional quality machine that provides a secure seal for EOGas sterilization bags. This rugged sealer features adjustable independent seal and cooling timers, and a built-in cutter knife.

Sealer Specifications

- ▶ Seal length: 24 1/2" Seal width: 1/8"
- ▶ Dimensions: 29"w x 6"h x 7"d
- ▶ 115 Volt standard, 220 Volt optional
- ▶ Power(watts): 575
- ▶ Shipping weight: 47lbs.
- ▶ 1 year warranty

Optional Sealer Stand is made of heavy gauge stainless steel, and is designed to support heavy EOGas bags without tipping.

Sealer Stand features a pull out shelf to hold EOGas bags during sealing.



AN1091
shown with
optional AN5099 stand

Stand Dimensions:

30.125"W x 28.125"D x 33.5"H
(76.51cm x 71.38cm x 85.1cm)

Extending shelf:

12"W x 16.5"L
(30.48cm x 41.91cm)

PN259_033009

Manufactured by: Andersen Sterilizers, Inc.



www.anpro.com
Andersen Products, Inc.
 Health Science Park • 3202 Caroline Drive
 Haw River, NC 27258-9564 USA
 800-523-1276 • fax 336-376-8153

www.andcal.com
H. W. Andersen Products of California, Inc.
 Health Science Park • 3151 Caroline Drive
 Haw River, NC 27258-9789 USA
 800-524-3455 • fax 336-376-3088

Distributed in Europe by:
H.W Andersen Products, Ltd.
 Davy Road • Clacton-on-Sea
 Essex CO15 4XA UK

Manufactured by:
Andersen Sterilizers, Inc.
 3154 Caroline Drive
 Haw River, NC 27258 USA